

NEWS

doctor smile[®]
dental laser



doctor smile goes WORLDWIDE!

November 18th, 2011

The word over doctor smile dental lasers is quickly spreading across the globe: all over the world wide web our Simpler laser is being noticed, from India, across the middle east to America! Find out more about the advantages of laser dentistry at www.doctor-smile.com



doctor smile[®]
dental laser

Home Leaders Buyers Guide News Articles Products Videos Online CE Meetings FREEInfo

Search...

Today's Dental News

> Today's Dental News

Industry News >

> Events

> Science and Medicine

Doctor Smile Introduces New Laser

Wednesday, 27 July 2011 14:45

Doctor Smile, one of the world's leading dental laser company, is pleased to announce the latest Diode Laser product: **Simpler laser**

The Doctor Smile Simpler 8 Watt diode laser with built-in NAVIGATION GUIDE which will carefully escort the doctor to perform a lot of operations with greater In Surgery, the T.O.P. - Tissue Optimized Pulsing configuration, the doctor can easily control the Simpler laser beam for the three soft tissue biotype: Fibrous, Normal and Granulation. TOP optimizes cutting and coagulation with a bloodless, painless and quicker healing operations. Simpler laser is an indispensable tool in the daily routine of dental practice. The laser can be used for an incredible number of applications that work alongside, and in many cases replace, traditional techniques and therapeutic protocols.

The Doctor Smile laser works in full respect of the tissue present in the oral cavity, therefore it is precise, safe, efficient, and in the majority of cases painless. Using a laser reduces the possible need of anaesthetic, sutures and pharmaceuticals when carrying out surgical and therapeutic interventions. A fast recovery therefore, which avoids unpleasant post-operative complications.

Simpler laser is an indispensable tool in the daily routine of dental practice.



DENTAL360 Dental Product Shopper Hygiene Product Shopper Dental Learning.net Forum From The Podium Practice Building Articles Blog About Us Contact Us Register

DENTAL 360

THEIR VOTE THEIR LASER



Dental PRODUCT SHOPPER YOUR SOURCE FOR PRODUCT EVALUATIONS & INFORMATION

DENTAL PRODUCTS

HYGIENE PRODUCTS

DENTAL LEAR

Home > Dental Products > Simpler Laser

Product Evaluations New Products Promotions Manufacturer Index Equ

Categories: Emerging Technologies in Dentistry | Lasers | Lasers

SIMPLER Laser

The Doctor Smile SIMPLER 8 W diode laser with built-in NAVIGATION GUIDE provides the clinician with efficient operation through a list of preset treatments with parameters and clinical information. In Surgery mode with the Tissue Optimized Pulsing (TOP) configuration, the dentist can easily control the SIMPLER laser beam for the 3 soft-tissue biotypes: fibrous, normal, and granulation.

TOP optimizes cutting and coagulation with a bloodless and quicker healing. The SIMPLER laser can be used for a great number of applications that work alongside, and in many cases replace, traditional techniques and therapeutic protocols. The Doctor Smile SIMPLER laser works in the tissue present in the oral cavity in a precise, safe, and efficient manner, and in the majority of cases, is also painless for the patient.

Request More Info



Dental NEWS

PRODUCT NEWS

Home Product News Magazine Calendar Events Links Media Kit Subscribe Contact Us



click

DOCTOR SMILE SIMPLER

The Doctor Smile Simpler 8 Watt diode laser works in full respect of the tissue in the oral cavity, therefore it is precise, safe, efficient, and in the majority of cases painless. Moreover the built-in NAVIGATION GUIDE will carefully escort the doctor to perform a lot of operations with greater efficiency. Using a laser reduces the possible need of anaesthetic, sutures and pharmaceuticals when carrying out surgical and therapeutic interventions. Laser treatment can be used to treat multiple pathologies such as gum care, desensitisation, surgery, whitening and aesthetics dentistry. In Surgery, the T.O.P. - Tissue Optimized Pulsing configuration, the doctor can easily control the Simpler laser beam for the three soft tissue biotype: Fibrous, Normal and Granulation. TOP optimizes cutting and coagulation with a bloodless, painless and quicker healing operations. Simpler laser is an indispensable tool in the daily routine of dental practice.

contact: info@doctor-smile.com

Details at: www.doctor-smile.com



Copyright © 2011 dentalnews.com